

Series OXY-SEN Percent Oxygen Monitor

- Economically priced
- Standard measuring range of 0-100%
- Easy to read display
- Two 4-20ma analog outputs scaled at 0-25% and 0-100%
- Open diffuser sensor designed to measure the contents of a chamber without a sampling system
- Easy to calibrate on room air

Trace PPM level alternatives also available



Key applications

Glovebox Chambers



Anaerobic Chambers



Hypoxia Chambers



PREMIUM INSIGHTS – GAIN REAL-TIME INSIGHT INTO YOUR PROCESS

Process Insights' products and solutions deliver innovative and differentiated analysis and measurement solutions and technologies that add high value to our customers and protect the environment.

Our commitment is to deliver smart and affordable innovation that optimizes process, improves safety, and transforms our world.

CENTERS OF EXCELLENCE

PROVIDING PROVEN SOLUTIONS FROM A GLOBAL TECHNOLOGY LEADER

COSA Xentaur

4140 World Houston Parkway
Suite 180
Houston, TX 77032
USA
+1 713 947 9591
www.alphaomegainstruments.com
sales@cosaxentaur.com

Process Insights – The Americas

4140 World Houston Parkway
Suite 180
Houston, TX 77032
USA
+1 713 947 9591
info@process-insights.com

Process Insights – EMEA

ATRICOM
Lyoner Strasse 15
60528 Frankfurt
Germany
+49 69 20436910
info@emea.process-insights.com

Process Insights – APAC

Wujiang Economic and Technology
Development Zone
No. 258 Yi He Road, 215200 Suzhou
Jiangsu Province China
+86 400 086 0106
info@apac.process-insights.com

OUR PREMIUM GLOBAL BRANDS

